

Skin Grafting Family

The Power of Precision



Our Range of Best Selling Instruments

Surgical Instruments ■ Chirurgische Instrumente ■ Instrumentos para quirúrgicos
Instruments de chirurgie ■ Strumenti chirurgici

 +92 300 0411477

 www.sterimed-instruments.com

 sterimedinstruments@gmail.com

 Harrar P.O Khas, Tehsil & District Sialkot - 51310 - Pakistan




Dynamic Meshing Options

SteriGraft™ II Tissue Expansion System

Using the grooved side of the Skin Graft Dermacarriers™, the SteriGraft II System offers four expansion ratios – 1.5:1, 3:1, 6:1, and 9:1.

The SteriGraft II instrument consists of a ratchet handle, continuous feed roller, a cutting roller, and a guidance plateau. The guidance plateau helps ensure proper alignment of the Dermacarrier II Skin Graft Carriers and cutting rollers. The ratchet handle facilitates proper advancement of the carrier.



SteriGraft II Tissue Expansion System

With a variety of ratios offered by the Sterimed Meshers, you can reduce clutter in the OR by minimizing the amount of equipment needed during a procedure.



Sterimed Skin Graft Mesher Carriers

Dermacarriers Length 8 in. (20.3 cm)

Sterimed II Tissue Expansion System

Sterigraft II Tissue Expansion System Complete	SMSG-200-00
Sterigraft II Mesher	SMSG-200-10
Ratchet Handle	SMSG-200-20

SteriGraft II Carriers

1.5:1 Dermacarriers II Skin Graft Carrier	SMSG-210-15
3:1 Dermacarriers II Skin Graft Carrier	SMSG-210-30
6:1 Dermacarriers II Skin Graft Carrier	SMSG-210-60
9:1 Dermacarriers II Skin Graft Carrier	SMSG-210-90

Service and Repair Programs

All Sterimed dermatomes are recommended for service and electronic calibration annually.




Packaging:

Non-Sterile Instruments in Poly Bag &

EO Sterile Instruments in Sterilized Pouch

3:1 Ratio
Non-Sterile



6:1 Ratio
Non-Sterile



ETO Packing

